

PROTEIN in CSF

Folin-Lowry Colorimetric Method 100 Tests

PRINCIPLE :

Two reactions are involved: (a) an initial interaction of protein and Cu^{2+} in alkali (related to biuret reaction); (b) a reduction of the phosphotungstic and phosphomolybdic acids to molybdenum blue and tungsten blue both by the Cu-protein complex and by the tyrosine and tryptophan of the protein. The latter two give color in the absence of Cu^{2+} , but the rest of the protein gives no color without Cu^{2+} . About 75% of the color is dependent on the Cu^{2+} . The maximum absorption of the colored products is at 750 nm.

SAMPLE :

CSF

Normal Values :

15 - 45 mg / dl

REAGENTS :

1-	Standard	50 mg / dL
2-	Alkaline Copper - Tartrate	
	Sodium Carbonate	200 mmol / L
	Sodium Tartrate	20 mmol / L
	Sodium Hydroxide	800 mmol / L
	Copper Sulphate	6 mmol / L
3-	Folin and Ciocalteu phenol Reagent	

STABILITY :

The reagents are stable up to the expiry date specified when stored at +4 to +8 °C away from light .

REFERENCE :

Daughaday WH, Lowry OH, Rosebrough NJ, Fields Ws. J, Lab. Clin. Med. 39, 663, (1952).

PROCEDURE:

	Blank (ml)	Standard (ml)	Sample (ml)
D. Water	0.05	-	-
Standard	-	0.05	-
Sample	-	-	0.05
Reagent 2	1.0	1.0	1.0
Mix well, let Stand 15 min. add rapidly:			
Reagent 3	0.05	0.05	0.05

Mix well immediately. Incubate for 10 min. at 37°C. Read absorbances at 700 (650 – 750 nm) for sample (A_{Sample}) and standard (A_{standard}) against reagent blank. Color stable for Several hours.

CALCULATION :

$$\text{Protein concentration (mg/dL)} = \frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times 50$$

PROTEIN in CSF

Folin – Lowry Colorimetric Method

+4 °C 100 Tests
In vitro diagnostic use

CAT. NO. TP 20 21

REAGENTS

R1	Standard	2	ml
R2	Alkaline reagent	100	ml
R3	Phenol reagent	5	ml

www.bio-diagnostic.com

29 Tahreer St., Dokki, Giza, Egypt

TeleFax : 3385184

e.mail : biodiagnostic_eka@lycos.com

info@bio-diagnostic.com

customerservice@bio-diagnostic.com



PROTEIN in Urine

Folin-Lowry Colorimetric Method 100 Tests

PRINCIPLE :

Two reactions are involved: (a) an initial interaction of protein and Cu^{2+} in alkali (related to biuret reaction); (b) a reduction of the phosphotungstic and phosphomolybdic acids to molybdenum blue and tungsten blue both by the Cu-protein complex and by the tyrosine and tryptophan of the protein. The latter two give color in the absence of Cu^{2+} , but the rest of the protein gives no color without Cu^{2+} . About 75% of the color is dependent on the Cu^{2+} . The maximum absorption of the colored products is at 750 nm.

SAMPLE :

Urine

Normal Values :

40 -400 mg/ 24 hrs

REAGENTS :

1-	Standard	50 mg / dL
2-	Alkaline Copper - Tartrate	
	Sodium Carbonate	200 mmol / L
	Sodium Tartrate	20 mmol / L
	Sodium Hydroxide	800 mmol / L
	Copper Sulphate	6 mmol / L
3-	Folin and Ciocalteu phenol Reagent	

Additional Reagents:

- Trichloro acetic acid (TCA) 20 g/dL
- Sodium hydroxide (NaOH) 1 g/dL

STABILITY :

The reagents are stable up to the expiry date specified when stored at +4 to +8 °C away from light .

REFERENCE :

Daughaday WH, Lowry OH, Rosebrough NJ, Fields Ws. J, Lab. Clin. Med. 39, 663, (1952).

PROCEDURE:

STEP FOR ELIMINATION THE INTERFERENCE:

Precipitation of protein:

Urine	2.0 ml
TCA 20%	0.5 ml

Mix well. Leave for 10 min. at R. Temp. Centrifuge for 15 min. at 3000 rpm. Discard supernatant. Dissolve ppt in 0.05 ml NaOH (1g/dL) then add 1.95 ml d. water, mix Well, use for assay:

	Blank (ml)	Standard (ml)	Sample (ml)
D. Water	0.05	-	-
Standard	-	0.05	-
Sample	-	-	0.05
Reagent 2	1.0	1.0	1.0
Mix well, let Stand 15 min. add:			
Reagent 3	0.05	0.05	0.05

Mix well immediately. Incubate for 10 min. at 37°C. Read absorbances at 700 (650 - 750 nm) for sample (A_{Sample}) and standard (A_{standard}) against reagent blank. Color stable for several hours.

CALCULATION :

$$\text{Protein Concentration (mg/dL)} = \frac{A_{\text{sample}}}{A_{\text{standard}}} \times 50$$

PROTEIN in Urine

Folin – Lowry Colorimetric Method

+4 °C 100 Tests
In vitro diagnostic use

CAT. NO. TP 20 21

REAGENTS

R1	Standard	2 ml
R2	Alkaline reagent	100 ml
R3	Phenol reagent	5 ml

www.bio-diagnostic.com

29 Tahreer St., Dokki, Giza, Egypt

TeleFax : 3385184

e.mail : biodiagnostic_eka@lycos.com

info@bio-diagnostic.com

customerservice@bio-diagnostic.com

