

## PROTEIN in CSF

### Folin-Lowry Colorimetric Method 100 Tests

#### PRINCIPLE :

Two reactions are involved: (a) an initial interaction of protein and  $\text{Cu}^{2+}$  in alkali (related to biuret reaction); (b) a reduction of the phosphotungstic and phosphomolybdic acids to molybdenum blue and tungsten blue both by the Cu-protein complex and by the tyrosine and tryptophan of the protein. The latter two give color in the absence of  $\text{Cu}^{2+}$ , but the rest of the protein gives no color without  $\text{Cu}^{2+}$ . About 75% of the color is dependent on the  $\text{Cu}^{2+}$ . The maximum absorption of the colored products is at 750 nm.

#### SAMPLE : CSF

**Normal Values :**  
15 - 45 mg / dl

#### REAGENTS :

1-	Standard	50 mg / dL
2-	Alkaline Copper - Tartrate	
	Sodium Carbonate	200 mmol / L
	Sodium Tartrate	20 mmol / L
	Sodium Hydroxide	800 mmol / L
	Copper Sulphate	6 mmol / L
3-	Folin and Ciocalteu phenol Reagent	

#### STABILITY :

The reagents are stable up to the expiry date specified when stored at +4 to +8 °C away from light .

#### PROCEDURE:

	Blank ( ml )	Standard ( ml )	Sample ( ml )
D. Water	0.05	-	-
Standard	-	0.05	-
Sample	-	-	0.05
Reagent 2	1.0	1.0	1.0
Mix well, let Stand 15 min. add rapidly:			
Reagent 3	0.05	0.05	0.05

Mix well immediately. Incubate for 10 min. at 37°C. Read absorbances at 700 (650 – 750 nm) for sample ( $A_{\text{Sample}}$ ) and standard ( $A_{\text{standard}}$ ) against reagent blank. Color stable for Several hours.

#### CALCULATION :

$$\text{Protein concentration (mg/dL)} = \frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times 50$$

#### REFERENCE :

Daughaday WH, Lowry OH, Rosebrough NJ, Fields Ws. J, Lab. Clin. Med. 39, 663, (1952).

## **PROTEIN in CSF**

Folin – Lowry Colorimetric Method  
+4 °C 100 Tests  
In vitro diagnostic use

**CAT. NO. TP 20 21**

## **REAGENTS**

<b>R1</b>	<b>Standard</b>	<b>2 ml</b>
<b>R2</b>	<b>Alkaline reagent</b>	<b>100 ml</b>
<b>R3</b>	<b>Phenol reagent</b>	<b>5 ml</b>

## **CONTACTS**

Tele: 02-33385184

Mobil: 0109 – 349 20 77

Fax : 02-33385184 (102)

e.maile : [info@bio-diagnostic.com](mailto:info@bio-diagnostic.com)

Website: [www.bio-diagnostic.com](http://www.bio-diagnostic.com)

Adress: 29 Tahreer St., Dokki, Giza, Egypt