


TRIGLYCERIDES

Enzymatic Colorimetric Method 50 Tests

PRINCIPLE :



SAMPLE :

Serum or plasma collected in heparin, sodium fluoride, oxalate, citrate or EDTA. Hemolysis will interfere .

NORMAL VALUES :

Women : 40 - 140 mg. / dL (0.46 – 1.60 mmol/L)
Men : 60 – 165 mg / dL (0.68 – 1.89 mmol/L)

REAGENTS :

1.	Standard	200 mg / dL (2.29 mmol/L)
2.	Buffer – Chromogen:	
	Buffer pH 7.5	100 mmol /L
	4 -Chlorophenol	3 mmol /L
3.	Enzymes :	
	Lipase	> 1000 U / L
	Glycerokinase	> 400 U / L
	Glycerol-3- phosphate oxidase	> 2000 U / L
	Peroxidase	> 200 U / L
	4- Aminoantipyrine	0.5 mmol / L
	ATP	0.5 mmol / L

STABILITY :

The reagents are stable up to the expiry date specified when stored at +4 to +8 °C .

PROCEDURE :

Working reagent : Mix equal volumes of reagent 2 and 3 immediately before the assay .

	Blank ml	Standard ml	Sample ml
Standard	-	0.02	-
Sample	-	-	0.02
Working reagent	1.0	1.0	1.0

Mix well. Incubate for 10 min. at 37°C. Measure the absorbances of the sample (A_{Sample}) and the standard (A_{Standard}) against blank, at 505 nm. (492 - 550) . The color intensity is stable for 30 min. Linearity up to 1000 mg / dL (11.4 mmol/L).

CALCULATION :

Triglycerides Concentration

$$= \frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times \text{Standard Conc.}$$

QUALITY CONTROL:

For accuracy and reproducibility control: Assayed Multi-Sera Normal and Elevated.

REFERENCE :

Fassati P ., Prencipe. L (1982) Clin. Chem., 28. 2077.

TRIGLYCERIDES

Enzymatic Colorimetric Method
+4 to +8°C 50 Tests
In vitro diagnostic use

CAT. No. TR 20 30

REAGENTS

R1 Standard	2 ml
R2 Buffer- Chromogen	25 ml
R3 Enzymes	25 ml

CONACTS

Tele: 02-33385184

Mobil: 0109 – 349 20 77

Fax : 02-33385184 (102)

e.maile : info@bio-diagnostic.com

Website: www.bio-diagnostic.com

Adress: 29 Tahreer St., Dokki, Giza, Egypt